


Biomerix 3D Scaffold™

The **Biomerix 3D Scaffold™** is designed to mimic the nature and function of the extracellular matrix (ECM) for stem cell research and therapeutic applications in regenerative medicine

- ❖ **Chemical Formulation**
 - Polycarbonate polyurethane-urea thermoset chemistry
- ❖ **Fully Reticulated 3D Morphology**
 - Interconnected 3D network of cells and pores
 - Open-cell, porous structure with void content of >90%
- ❖ **Demonstrated Biocompatibility with Multiple FDA Clearances**
 - Full panel of ISO 10993 biocompatibility testing completed
 - Ten regulatory clearances in the US, EU, and Canada for soft tissue repair, orthopedic, and vascular applications
- ❖ **Ideal Platform for Stem Cell Therapy & Regenerative Medicine**
 - Applications in drug discovery, tissue engineering, & cell banking
- ❖ **Available Form Factors & Pricing**
 - 25 cylindrical discs, 5 mm diameter x 2 mm thickness
 - Other custom form factors available upon request



Biomerix 3D Scaffold™ Features & Benefits

Features	Mimic <i>In-Vivo</i> Conditions	Enhance Cell Viability	Promote ECM Synthesis	Optimize 3D Cell-ECM Interactions	Provide Stable Environment	Support Cell Differentiation	Maximize Throughput
Biostable Polyurethane Chemistry					✓		✓
Interconnected 3D Network of Cells & Pores	✓	✓	✓	✓		✓	✓
High (>90%) Porosity	✓	✓	✓	✓		✓	✓
High Surface Area to Volume Ratio	✓	✓	✓	✓		✓	✓
High Affinity for Proteins & Peptides	✓	✓		✓		✓	✓
Reproducible-High Volume Manufacturing Process		✓		✓	✓	✓	✓
Biocompatible	✓	✓	✓	✓	✓	✓	✓

Studies Evaluating Biomerix 3D Scaffold™ in Regenerative Medicine

Ex-Vivo Hematopoietic Stem Cell Niche in the Rotary Cell Culture System™

Collaboration with Dr. Stephen Navran, Synthecon Inc., TX

Objective:

- Create a bioartificial hematopoietic stem cell niche to expand umbilical cord blood stem cells

Methods:


- Bone marrow stromal cells (MSC) expressing GFP cultured on *Biomerix 3D Scaffold* in Rotary Cell Culture System
- Umbilical Cord Blood (UCB) cells labeled with DiIC₁₈-DS introduced into the bioreactor in suspension

Results:

- UCB cells homed to the stroma/scaffold bioartificial niche
- The niche construct was maintained for at least 2 weeks

Conclusions:

- The combination of *Biomerix 3D Scaffold* and a dynamic rotary cell culture bioreactor can create a bioartificial niche for the expansion of rare UCB stem cells



Bone Marrow MSC and Umbilical Cord Blood Cells

Bone Marrow Hematopoietic Niche Formation in an *In Vivo* Transgenic Mouse Model

Collaboration with Dr. Daniel Kraft, Stanford University School of Medicine, CA

Objective:

- Investigate bone marrow hematopoietic niche formation in the *Biomerix 3D Scaffold*

Methods:


- Marrow cells isolated from adult and fetal mice (Luciferase transgenic)
- Artificial marrow created in the *Biomerix 3D Scaffold*
- Controls fabricated using intact whole bone marrow and Matrigel
- Constructs implanted into the renal capsule of GFP+ transgenic mice
- Donor cells tracked *in vivo* for up to eight months; homing of host cells to the niche tracked using GFP+

Results:

- Excellent results in terms of *in-vivo* viability of HSCs in the *Biomerix Scaffold*, with kidney capsule *Biomerix* grafts maintained *in-vivo* for up to eight months
- *Biomerix 3D Scaffold* allows transplanted HSCs to establish and expand
- *Biomerix 3D Scaffold* also formed osteogenic niche from the fetal cells

Conclusions:

- *Biomerix 3D Scaffold* is an excellent platform for bone marrow niche formation and expansion of hematopoietic stem cells



Kidney with *Biomerix 3D Scaffold* two weeks after implantation. Note strong bioluminescent signal from donor derived marrow cells. *Biomerix 3D Scaffold* appeared to be somewhat hard, consistent with early bone formation (osteogenesis). Cells from the scaffold were transplanted to secondary recipient after sublethal irradiation (400Gy). Engraftment is evident (middle mouse), from *Biomerix 3D Scaffold* derived transplanted cells.

Ectopic Human Mesenchymal Stem Cell-Coated Scaffolds in NOD/SCID Mice: An *In-Vivo* Model of the Leukemia Niche

Collaboration with Dr. Sarah R. Vaiselbuh, Feinstein Institute for Medical Research, NY

Objective:

- Evaluate potential of *Biomerix 3D Scaffold* to support malignant hematopoiesis within the human stromal microenvironment

Methods:


- MSCs isolated from human bone marrow obtained electively after orthopedic surgery
- MSC-coated *Biomerix 3D Scaffolds* prepared *in-vitro* for 24-48 hours
- Ectopic MSC-coated scaffolds implanted in NOD/SCID mice and injected with primary AML cells
- Known leukemia CXCR4 antagonist used to treat AML niche *in-vivo*
- Scaffolds harvested at 1, 4, 8, and 16 weeks and analyzed histologically

Results:

- Leukemia niche successfully maintained for 4 months in murine model, supporting a human bone marrow stromal network of human AML cells
- CXCR4 leukemia antagonist shown to disrupt MSC stromal formation

Conclusions:

- *Biomerix 3D Scaffold* provides excellent platform to create *in-vivo* disease models that can be used for targeted drug discovery



Novel Flow-Perfusion Bioreactor with *Biomerix 3D Scaffold*™

Collaboration with Dr. Stephen W. Warren, New York University Medical Center, NY

Objective:

- Demonstrate that large 3D constructs of the *Biomerix 3D Scaffold* support viable cellular distribution, density, and metabolic activity using murine preosteoblasts in a novel flow-perfusion bioreactor

Methods:

- Murine preosteoblastic cells seeded onto large *Biomerix 3D Scaffold* discs (24 x 6 mm)
- Static culture scaffolds maintained in culture plates
- Dynamic culture scaffolds loaded into custom flow-perfusion bioreactor designed to optimize chemotransportation
- Scaffolds harvested at 0, 2, 4, 6, and 8 days and analyzed histologically

Results:

- Flow-perfusion culture promotes cellular viability in the *Biomerix 3D Scaffold* with significantly higher cellular density and metabolic activity compared to static culture

Conclusions:

- Large constructs of the *Biomerix 3D Scaffold* support chemotransportation & viable cellular activity in a flow-perfusion bioreactor, and may enable prefabrication of biological constructs large enough to solve clinical problems

